

Dunaliella Cultivation

A Abubakar, R Swamy, P Harvey, University of Greenwich

2012



Grant Contract: FED/2009/217066

Determination of glycerol


There are two methods to measure glycerol, a chemical and an enzymatic method.

Chemical Method

Chen et al., (2011).

- Harvest the cells by centrifugation at 5,000 g for 15 min at room temperature.
- Decant the supernatant and re-suspend the algal pellet by centrifugation in 1.0 mL of distilled water and 0.2 mL of chloroform.
- Sonicate the suspension for 8 minutes in an ultrasonic cell disruptor at room temperature.
- Remove cell debris by centrifuging the suspension at 10,000g for 10 min at room temperature and collect the supernatant.
- Take 100 µl of supernatant and make up to 0.2 mL with distilled water in a test tube.
- Add 1 mL of sodium periodate reagent (3 mmol/L sodium periodate and 100 mmol/L ammonium acetate in 100 mL 6% acetic acid) with mixing.
- After 5 minutes add 2.5 mL of acetyl acetone reagent (acetyl acetone/isopropanol=1:99) with mixing.
- Place the test tube in a water bath at 60 °C for 30 min. Place a marble at the mouth of the tube to restrict evaporative loss during incubation.
- Cool to room temperature then read the absorbance at 410 nm. The corresponding concentrations of blank samples without algae should be run in each test.

The absorbance of each sample can be converted to glycerol content by reference to a standard curve for pure glycerol (0-50 µg glycerol /0.2 ml) .



Dunaliella Cultivation

A Abubakar, R Swamy, P Harvey, University of Greenwich

2012



Grant Contract: FED/2009/217066

Enzymatic method


Microalgal extract preparation

After Ben-Amotz and Avron (1973).

- Harvest the algal cells by centrifugation at 5,000 g for 15 min at room temperature.
- Decant the supernatant and re-suspend in 1.5 mL of distilled water.
- Sonicate the suspension for 8 minutes in an ultrasonic cell disruptor at room temperature.
- Remove cell debris by centrifuging the suspension at 10,000g -12000g for 10 min at room temperature and collect the supernatant.
- Deproteinize the suspension by adding trichloroacetic acid to 2%w/v final concentration, then neutralize with NaHCO.

Enzymic Determination of glycerol

The glycerol content in the algal cells is determined by enzyme-coupled assay in which glycerol is oxidised to dihydroxyacetone and NAD⁺ reduced to NADH using glycerol dehydrogenase in the following reaction:



One unit of glycerol dehydrogenase reduces one micromole of NAD per minute at 25°C and pH 10.0 under the specified conditions.

Pipette into a 1ml UV cuvette the following:

- algal extract
- 0.1M NAD
- make up to 1 ml with 0.125 M carbonate/bicarbonate buffer, pH 10.0.

Start the reaction with the addition of 10 µl of glycerol dehydrogenase (100unit/ml in 0.05 M potassium phosphate buffer solution).

Determine the increase in absorbance at 340 nm resulting from the reduction of NAD. The concentration of glycerol can be calculated using the extinction coefficient for NADH ($\epsilon_{340} = 6220 \text{ M}^{-1}\text{cm}^{-1}$)

References

1. Johnson, M.K., Johnson, E.J, Mac Elroy, R.D., Speer, H.L. and Bruff, B.S. (1968) Effects of salts on the halophilic alga *Dunaliella viridis*, J. Bacteriology 95: 1461-1468.
2. Borowitzka MA (1988) Vitamins and fine chemicals. In Borowitzka MA, Borowitzka LJ (eds), Microalgal Biotechnology. Cambridge University Press, Cambridge, 153-196.

Dunaliella Cultivation

A Abubakar, R Swamy, P Harvey, University of Greenwich

2012



Grant Contract: FED/2009/217066

3. Chen, H., Lao, Y.M., Jiang, J.G. (2011). Effects of salinities on the gene expression of a (NAD⁺)-dependent glycerol-3-phosphate dehydrogenase in *Dunaliella salina*. *Sci.Total Environ.*, 409:1291–7

Other Useful References

1. Ben-Amotz A. & Avron M. 1973 The role of glycerol in the osmotic regulation of the halophilic alga *Dunaliella parva*. *Plant Physiol.* **51**, 875-878 (DOI: 10.1104/pp.51.5.875)
2. Avron M., Ben Amotz A. 1978 Production of glycerol from algae. *US Patent no. 4115949*
3. Avron M., Ben Amotz A. 1980 Production of glycerol, carotenes and algae meal. *US Patent no. 4199895*
4. Ben-Amotz A., Shaish M. & Avron M. 1991. The Biotechnology of Cultivating *Dunaliella* for Production of β -Carotene Rich Algae. *Biosource Technol* **38**, 233-235 (DOI:10.1016/0960-8524(91)90160-L)
4. Ben-Amotz A., Polle E.W. & Subba Rao D.V. 2009. (Eds) *The Alga Dunaliella: Biodiversity, Physiology, Genomics and Biotechnology*. Science Publishers, Jersey, Plymouth, 555pp
5. Ben-Amotz A., Avron M. 1990. The biotechnology of cultivating the halotolerant alga *Dunaliella*. *Tibtech* **8**, 121-126 (DOI:[10.1016/0167-7799\(90\)90152-N](https://doi.org/10.1016/0167-7799(90)90152-N))
6. L. J. Borowitzka, M. A. Borowitzka, T. P. Moulton (1984) *Hydrobiologia* 116-117, pp 115-121 The mass culture of *Dunaliella salina* for fine chemicals: From laboratory to pilot plant